
INFORMATION LEAFLET
NOTICE EXPLICATIVE
PACKUNGSBEILAGE
INDLÆGSSEDDEL
BIJSLUITER
FOLLETO ILLUSTRATIVO
FOGLIO ILLUSTRATIVO
INFORMATIONSBLAD FOLHETO
ILLUSTRATIVO ΦΥΛΛΑΔΙΟ
ΕΠΕΞΗΓΗΣΕΩΝ
説明書

- * Refers to IIa class devices
- * Concerne les dispositifs de classe IIa
- * Es bezieht sich auf Geräte der Klasse IIa
- * Henvist til aparater klasse IIa
- * Betreft klasse IIa
- * Se refiere a dispositivos de clase IIa
- * Riferito ai dispositivi classe IIa
- * Avsågse device klass IIa
- * Refere-se a dispositivos de classe IIa
- * Αναφέρεται σε συσκευές κατηγορίας IIa
- * デバイスクラスIIa

English:	5	EN
Français:	13	FR
Deutsch:	21	DE
Dansk:	29	DA
Nederlands:	37	NL
Español:	45	ES
Italiano:	53	IT
Svenska:	61	SE
Português:	69	PT
Ελληνικά:	77	GR
日本語:	85	JP

ENGLISH
INFORMATION LEAFLET
HPF S.r.l. INSTRUMENTS

WARNINGS

- Critical elements, such as handles, articulated joint instruments or threads and in particular all those medical instruments with atraumatic toothing require particularly careful inspections.
- In order to avoid damage and subsequent corrosion due to rubbing of metals, avoid under all circumstances the use of metal brushes or metal sponges for removing stains.
- Do not exceed 137 °C.
- Devices that cannot be disassembled or that are complex require particular attention during cleaning.
- After having disinfected and cleaned the instruments, immediately dry them.
- In case of thermolabile material, steam sterilization is preferable.
- Machine treatment with thermal disinfection and steam sterilization should be applied.
- Acetabular cutters: do not use hand-pieces with a torque higher than 25 Nm and the axial thrust is allowed to a single user.
- Flexible mandrel: do not subject the instrument to a torque higher than 7 Nm.

DESCRIPTION

HPF S.r.l. instruments are medical instruments used for surgeries:

- manual medical instruments belonging to class I;

- energy transfer instruments belonging to class II a (Refer to the Directive 93/42/EEC). **HPF** S.r.l. medical devices are authorized for use only by specialized, trained personnel certified to perform surgeries or handle related surgical instruments.

USE

Before clinical use, the surgeon must be familiar with all aspects of the surgery and the limits of the instruments.

These instruments have a duration time limit; this limit is determined by the level of wear and tear and damage arising from the intended use.

Therefore, it is recommended to replace the instruments when the user notices reduced performance or in the case visible signs of wear and tear appear.

RECOMMENDATIONS FOR USE

Generally, device treatment includes:

- Preparation (pre-treatment, gathering of instruments, pre-cleaning and disassembly, if applicable).
- Cleaning, disinfection, rinsing and drying.
- Visual inspection of cleaning and perfect conditions of the material.
- Verification of correct operation.
- Packaging and sterilization, authorization and storing, whenever applicable. After use, the instruments must be properly conditioned.

MANUAL CLEANING

- Reusable instruments must be cleaned immediately after use.
- Recesses (cavities) and hidden areas within an instrument must be inspected periodically in order to ensure the absence of residual substances.
- For manual cleaning, we recommend using a soft cloth with low particle release or synthetic brushes.
- After manual disinfectant cleaning it is necessary to perform intense and thorough rinsing using clear running water. During the procedure, manually remove any possible dirt residue. Then proceed to a final rinse using demineralized water and thoroughly dry the instruments. A good cleaning of the instruments plays also an important role in maintaining the value of the instruments and sets the stage for a successful sterilization.
- Pursuant to the international standards (EN ISO 15883) and national guidelines, only validated procedures for cleaning and disinfecting should be applied.

MECHANICAL CLEANING

- Reusable instruments must be cleaned immediately after use.
- Recesses (cavities) and hidden areas within an instrument must be inspected periodically in order to ensure the absence of residual substances.
- The best method to achieve standardization, cleaning and

disinfection is to use the machine treatment method. A good cleaning of the instruments plays also an important role in maintaining the value of the instruments and sets the stage for a successful sterilization.

- Pursuant to the international standards (EN ISO 15883) and national guidelines, only validated procedures for cleaning and disinfecting should be applied.

STERILIZATION

- This instrument exits the premises of **HPF** S.r.l. non-sterilized.
- All instruments and parts of instruments must be sterilized before use.
- The use of sterilized instruments on the patient establishes, within the framework of the European regulations in force that all the instruments are adequately cleaned and disinfected, sterilized in a suitable sterile packaging using a validated procedure and stored according to the regulations in force for sterile materials. This is why it is important to use only methods and sterilization devices that allow a validated sterilization process.
- All instruments must be removed from their packaging prior to sterilization.
- Reusable instruments must be disinfected and sterilized after use.
- The recommended sterilization for the instruments manufactured by **HPF** S.r.l. is steam sterilization (autoclave)

according to the EN ISO 17665 Directive and European pharmacopoeia regulations:

Sterilize by steam autoclaving, utilizing a pre-vacuum cycle **HPF** recommends the following cycle: Steam autoclave 132-135°C (270-275°F), minimum holding time 10 minutes.

- Any other validated pre- vacuum autoclave cycle may be used in alternative.
- Processed items should be stored in clean and secure stores to avoid damage or tampering.
- **HPF** S.r.l. accepts no responsibility in case problems arising from sterilizations different from the recommended one.
- The instruments must not come into contact with fluoridated or chloric treated products or with greasy detergents.
- Instruments with synthetic material elements (plastics) must never come in contact with strong acids (pH'4), with bases, organic solvents or ammonia, and not even with oxidized chemical products or any type product that may alter the materials.

PROHIBITION

It is strictly forbidden to apply any modification to the instruments manufactured by **HPF** S.r.l.; **HPF** S.r.l. is the only party with the authority to make such modifications. In case of failure to heed this prohibition, **HPF** S.r.l. accepts no responsibility for the consequences.

STORAGE

EN

- The instruments are sensitive to scratches and denting. Scratches on their surface can increase the risk of corrosion. For this reason it is recommended to handle them with care.
- Storage of medical instruments must be in an environment free from moisture in order to avoid excessive corrosion. This is valid also for the transportation and conditioning of the instruments.
- The surgical instruments are subject to wear and tear; thus, prior to use it is necessary to inspect them to verify that they are free from defects.

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ES

The first part of the paper discusses the importance of the research and the objectives of the study. It then presents a literature review of the existing research on the topic. The second part of the paper describes the methodology used in the study, including the data collection and analysis techniques. The third part of the paper presents the results of the study, and the fourth part discusses the implications of the findings. The paper concludes with a summary of the main findings and a list of references.

The research was conducted in a systematic and rigorous manner, following the principles of good research practice. The data was collected from a representative sample of the population, and the analysis was carried out using appropriate statistical methods. The results of the study are presented in a clear and concise manner, and the implications of the findings are discussed in detail. The paper is well-written and easy to read, and it provides a valuable contribution to the field of research.

The findings of the study have important implications for the field of research, and they provide a basis for further research. The research has shown that there is a need for more research in this area, and it has identified some of the key areas for future research. The research has also shown that there are some important factors that influence the results of the study, and it has provided some suggestions for how these factors can be controlled in future research.

The research has been carried out in a transparent and open manner, and the results of the study are available to the public. The research has been peer-reviewed, and it has been published in a reputable journal. The research is a valuable contribution to the field of research, and it provides a basis for further research.

The research has been carried out in a systematic and rigorous manner, following the principles of good research practice. The data was collected from a representative sample of the population, and the analysis was carried out using appropriate statistical methods. The results of the study are presented in a clear and concise manner, and the implications of the findings are discussed in detail. The paper is well-written and easy to read, and it provides a valuable contribution to the field of research.

The findings of the study have important implications for the field of research, and they provide a basis for further research. The research has shown that there is a need for more research in this area, and it has identified some of the key areas for future research. The research has also shown that there are some important factors that influence the results of the study, and it has provided some suggestions for how these factors can be controlled in future research.

The research has been carried out in a transparent and open manner, and the results of the study are available to the public. The research has been peer-reviewed, and it has been published in a reputable journal. The research is a valuable contribution to the field of research, and it provides a basis for further research.

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